Technovit® 7200/7210/7230 VLC

D Gebrauchsinformationen Technovit 7200/7210/7230 VLC

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GB User instructions Technovit 7200/7210/7230 VLC

Technovit 7200 VLC
Photocuring one-component methacrylate-based resin for embedding and subsequent production of thin ground sections in medicine and dentistry.

Technovit 7210 VLC
Photocuring one-component precision adhesive on the basis of methacrylate for precision adhesion of polymerized specimens on acrylic glass specimen holders.

Technovit 7230 VLC
Photocuring one-component adhesive methacrylate-based, for the fixing of tissue specimens in the embedding mould.

Properties
Technovit 7200 VLC was specially developed for the thin grinding technique. It completely penetrates hard tissue, and with this technique no decalcifi-cation of the tissue is necessary. Resins like compo-site tooth filling materials or bone cements are not affected by Technovit 7200 VLC.
Polymerization of the embedding medium takes place in the Histolux light unit. When adhering to the irradiation times indicated, temperatures will not exceed 40 ° C.
Its mechanical strength makes Technovit 7200 VLC particularly suitable for the thin grinding technique, and the sawing and grinding equipment manufac-tured by Exakt is best suited for this method.
All customary dyeing methods may be applied with high-quality results and in a minimum of time.

Handling
Fixation
Any fixation method customary in light microscopy may be used for this photocuring embedding sys-tem.

Dehydration
Dehydration of the tissue may be done via an ascending glycolmethacrylate or alcohol range.

Example
Glycolmethacrylate/Water = 70/30 3 hours
Glycolmethacrylate/Water = 80/20 3 hours
Glycolmethacrylate/Water = 90/10 3 hours
Glycolmethacrylate =100 3 hours
Glycolmethacrylate =100 3 hours
Glycolmethacrylate =100 3 hours

Infiltration
Infiltration is initiated with a mixture of alcohol/ Technovit 7200 VLC = 1/1 in the alcohol range or glycolmethacrylate/Technovit 7200 VLC = 1/1 in the glycolmethacrylate range, followed by infiltration with pure Technovit 7200 VLC.
If dehydration is done mechanically and the infiltra-tion with agitation and vacuum, dehydration times of 3 hours per stage and infiltration times of 8 hours per stage will suffice for tissues of 2 to 3 mm thick-ness. For dehydration and infiltration without agita-tion and vacuum the times per stage should be mul-tiplied by four.
Infiltration should take place in light-proof containers or equipment as contact with light will start the poly-merization.

Embedding
For planeparallel embedding the prepared and infil-trated tissue specimen is placed in the embedding mould I or II on a drop of fixation adhesive Technovit

7230 VLC and gently pressed down. It is important that the surface to be examined should be on the bottom of the mould. Then Technovit 7200 VLC is poured onto the specimen in such a manner that the specimen does not float and no air bubbles are entrapped.

Polymerization
The polymerization of the embedding takes place in the photopolymerization unit Histolux in two stages. During the first stage the embedding medium is po-lymerized almost completely at a low light intensity so that the polymerization temperature does not ex-ceed 40° C and no cracks will form. The final poly-merization of the embedding material infiltrated in the tissue is done at high light intensity during the second stage. Total polymerization time is 6 hours max.

If a maximum of 40 ml embedding medium is used at an ambient temperature of 23 ° C, polymerization temperature will not exceed 40 ° C.

Precision adhesion
The precision adhesion of tissue specimens to an acrylic glass specimen holder is done with the pho-tocuring precision adhesive Technovit 7210 VLC. The photopolymerization of the adhesive takes place in the precision adhesion press manufactured by Exakt and is finalized after 15 minutes irradiation time.

Dyeing
Possible with all customary dyeing methods.

Risks and Safety precautions
Product contains: 2-hydroxyethyl methacrylate, methyl methacrylate
R_11 Highly flammable.
R_36/37/38 Irritating to eyes, respiratory system and skin.
R_43 May cause sensitisation by skin contact.
R_51/53 Toxic to aquatic organisms, may cause long-term adverse effects in the aqua-tic environment.
S_9 Keep container in a well-ventilated place.
S_16 Keep away from sources of ignition – No smoking.
S_24 Avoid contact with skin.
S_26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
S_28 After contact with skin, wash immedi-ately with plenty of soap and water.
S_33 Take precautionary measures against static discharges.
S_37 Wear suitable gloves.
S_43 In case of fire, use sand, carbon dioxide or powdered extinguishing agent. Never use water.
S_57 Use appropriate container to avoid environmental contamination.

Storage advice:
Keep vessels always closed. Not store above room temperature. After expiry date the materials should not be used anymore. Keep away from reach of children.

Delivery units
Technovit 7200 VLC
1000 ml liquid
Technovit 7210 VLC
2 x 15 ml liquid
Technovit 7230 VLC
100 ml sirup

Our technical advice, whether verbal, in writing or by way of trials, is given in good faith but without war-ranty, and this also applies where proprietary rights of third parties are involved. It does not release you from the obligation to test the products supplied by us as and to their suitability for the intended proces-ses and uses. The application, use and processing of the products are beyond our control, and there-fore, entirely your own responsibility. Should, in spite of this, liability be established for any damage, it will be limited to the value of the goods delivered by us and used by you. We will, of course, provide con-sistent quality of our products within the scope of our General Conditions of Sale and Delivery.

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Heraeus Kulzer GmbH